

SAFE+ EXTRA LONG. EXTRA SAFE.

ADDED SPLASH PROTECTION – due to long cuff
BARRIER PROTECTION* – good Chemical Resistance
CONFIDENT GRIP – thanks to Finger Texturing

30 cm length provides added splash protection
Protects against certain chemicals*
Finger texturing offers a confident grip
Nitrile material is suitable for sensitive skin
Striking blue colour


EXTRA LONG. EXTRA SAFE.

The Sempercure® Safe+ effectively safeguards your hands when added barrier as well as splash protection are required. Thanks to its 30cm-long cuff, wrists and forearms are better protected, which makes it a suitable choice for glove users who work with risk patients or handle substances that can be harmful to the skin.

- **ENHANCED SAFETY:**
additional splash protection thanks to the extra-long cuff
- **CHEMICAL RESISTANCE:**
barrier protection against certain liquids and chemicals
- **WORK WITH CONFIDENCE:**
finger texturing and good tactile sensitivity
- **SUITABLE FOR SENSITIVE SKIN:**
nitrile material helps to prevent type-1 allergies

TECHNICAL DATA

TYPE	powder-free disposable examination glove		
MATERIAL	Nitrile Butadiene Rubber (NBR)		
COLOUR	blue		
GLOVE SHAPE	flat with rolled rim		
SIZE/OVERALL LENGTH as per EN 455-2	XS, S, M, L, XL / median 300 mm		
WALL THICKNESS measured in double layer in the palm area	0.20 mm		
BARRIER PERFORMANCE Tightness as per EN 455-1	AQL 1.5		
FORCE AT BREAK as per EN 455-2	≥ 8 N		
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years		
PACKAGING	in dispenser box: XS, S, M, L 100 gloves loose, XL 90 pieces in transport box: XS, S, M, L 1000 gloves loose, XL 900 pieces		
LABELLING	EN 455-1/-2/-3/-4, EN 420, EN 421, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5, ISO 15223-1		
PURPOSE as set out in MD Regulation (EU) 2017/745 as set out in PPE Regulation (EU) 2016/425	Medical device class I Single-use protective glove – PPE category III (protection against chemical substances for limited time)*		
ARTICLE NUMBERS	Size XS 8267 58031	Size S 8267 58033	
	Size M 8267 58035	Size L 8267 58037	
	Size XL 8267 58039		



CONTACT US!

Semperit Investments Asia Pte. Ltd.
8 Jurong Town Hall Road, #29-03/06 The JTC Summit,
Singapore 609434, Singapore
Tel.: +65 6274 4861 · Fax.: +65 6274 6977

Sempermed Europe GmbH
Am Belvedere 10,
A-1100 Vienna, Austria
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630

sempermed@semperitgroup.com · www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on <https://www.sempermed.com/knowledge-area/chemical-resistance/chemical-resistance-sempermed-brands/>
IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2022 Sempermed