

VELVET

SOFT AS VELVET AND LIGHT AS A FEATHER

YOUR DAILY CHOICE – a thin & light glove with a force at break of $\geq 6\text{ N}^*$
COMFORTABLE FIT – thanks to the soft glove design
ENHANCED SENSE OF TOUCH – good finger dexterity

High wearing comfort
Soft glove material

Enhanced sense
of touch



SEMPERCARE® VELVET

SOFT AND LIGHT AS A FEATHER

The comfortable and light Sempercure® Velvet is a soft latex-free option for professionals looking for high wearing comfort and an enhanced sense of touch combined with effective protection against infections.

- Feather light and velvety soft, the Sempercure® Velvet offers good wearing comfort and tactile sensitivity.
- The classification as medical device (class I) as set out in MD Regulation (EU) 2017/745 and as protective glove (category III) as set out in PPE Regulation (EU) 2016/425 makes the Sempercure® Velvet suitable for a variety of applications.
- Sempercure® Velvet is powder free (as per EN 455-3) and does not contain any natural rubber protein. As such, this glove is perfectly suitable for glove users with Type-I allergy.
- Individual size color coding on dispenser boxes allows for the easy identification of the correct glove size.
- Packed in 100 or 200 pieces per dispenser box.



TECHNICAL DATA

TYPE	unsterile, powder free examination glove for single use						
MATERIAL	nitrile Butadiene Rubber (NBR)						
COLOUR	lavender blue						
GLOVE SHAPE	Flat with rolled rim						
SIZE/OVERALL LENGTH as per EN 455-2	XS, S, M, L, XL / median 240 mm						
WALL THICKNESS measured in double layer in the palm area	min. 0.10 mm typical 0.11-0.12 mm						
BARRIER PERFORMANCE as per EN 455-1	AQL 1.5						
FORCE AT BREAK as per EN 455-2	median 6 N						
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years						
PACKAGING	available in 100 and 200 pieces per dispenser box in dispenser box XS, S, M, L 100 / XL 90 pieces in dispenser box XS, S, M, L 200 / XL 180 pieces 10 dispenser boxes per shipping carton						
LABELLING	EN 420, EN ISO 374-1, EN 374-2, EN 16523-1, EN 374-4, EN ISO 374-5, EN 455 1-4, EN ISO 15223-1						
PURPOSE as set out in MD Regulation (EU) 2017/745 as set out in PPE Regulation (EU) 2016/425	medical device CE class I single-use protective glove – PPE CE category III (protection against chemical substances for limited time)*						
ARTICLE NUMBERS	100 pcs.	size XS	300 000 6795	size M	300 000 6797	size XL	300 000 6799
		size S	300 000 6796	size L	300 000 6798		
	200 pcs.	size XS	300 000 6800	size M	300 000 6802	size XL	300 000 6804
		size S	300 000 6801	size L	300 000 6803		



CONTACT US!

Semperit Investments Asia Pte. Ltd.
8 Jurong Town Hall Road, #29-03/06 The JTC Summit,
Singapore 609434, Singapore
Tel.: +65 6274 4861 · Fax.: +65 6274 6977

Sempermed Europe GmbH
Am Belvedere 10,
A-1100 Vienna, Austria
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630

sempermed@semperitgroup.com · www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on <https://www.sempermed.com/knowledge-area/chemical-resistance/chemical-resistance-sempermed-brands/>
IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2022 Sempermed