

sempermed® 
supreme plus

SUPREME PLUS MORE GRIP TO YOUR WORLD

SECURE GRIP – confident instrument handling
SKIN-FRIENDLY – the Semperit Post Production Process
OPTIMUM SAFETY – ideal tactile properties



100 
+ YEARS
- EST. 1920 -
sempermed®
YOUR GLOBAL GLOVE EXPERT™

MORE GRIP TO YOUR WORLD

The Sempermed® Supreme Plus provides you with more grip in every situation while combining optimised safety, tactical properties and instrument guidance.

■ SKIN-FRIENDLY

The specially developed, automatic finishing process contains intensive and multiple washing steps. Skin-irritating substances resulting from the production process and the natural latex protein concentration are thus reduced to the technically lowest level.

■ THE LINING SYSTEM

The special net-like lining structure optimises the contact surface between the film and the skin, thus considerably reducing the friction resistance when putting on the glove. The result: The Sempermed® Supreme Plus facilitates easy donning of the gloves.

■ DOUBLE DONNING

For your safety and easy donning use Sempermed® Supreme as under glove and Sempermed® Supreme Plus as outer glove. Recommendation: If you wear different sizes, put on the larger pair first. Or use gloves of the same size.



TECHNICAL DATA

sempermed® supreme plus	
TYPE	Sterile surgical glove for single use, powder free, with synthetic lining
MATERIAL	Natural rubber latex (NR)
COLOUR	Natural white
GLOVE SHAPE	Fully anatomical with rolled rim
SIZE/OVERALL LENGTH as per EN 455-2	5½, 6 and 6½ min. 270 mm 7, 7½ and 8 min. 280 mm 8½ and 9 min. 285 mm
WALL THICKNESS measured in single layer in the palm area	0.19–0.23 mm
BARRIER PERFORMANCE Tightness exceeding EN 455-1 (AQL 0.65 instead of AQL 1.5)	AQL 0.65
FORCE AT BREAK exceeding EN 455-2 (≥ 12 N instead 9 N)	≥ 12 N
DURABILITY	3 years in original packaging if stored as per DIN 7716, ISO 2230
STERILISATION	Radiation STERILE R
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223
PURPOSE as set out in MD Directive 93/42/EEC as set out in PPE Regulation (EU) 2016/425	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for <u>limited time</u>)*
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical fibre-free peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs
ARTICLE NUMBERS	Size 5½ 822851521 Size 7½ 822851721 Size 6 822851601 Size 8 822851801 Size 6½ 822851621 Size 8½ 822851821 Size 7 822851701 Size 9 822851901



CONTACT US!

Semperit Investments Asia Pte. Ltd.
8 Jurong Town Hall Road, #29–03/06 The JTC Summit,
Singapore 609434, Singapore
Tel.: +65 6274 4861 · Fax.: +65 6274 6977

Sempermed Europe GmbH
Am Belvedere 10,
A-1100 Vienna, Austria
Tel.: +43 1 79 777-0 · Fax: +43 1 79 777-630

sempermed@semperitgroup.com · www.sempermed.com

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/chemical-resistance-brands

IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the «Products» tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. **Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used.** In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. **CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions.** ©Copyright 2022 Sempermed