

 **semper
med**

supreme plus

More grip to your world

The sempermed® supreme plus provides you with more grip in every situation while combining optimized safety, tactical properties and instrument guidance.



Easy to don



Skin-friendly



Secure grip



supreme plus



Easy to don

The special net-like lining structure optimizes the contact surface between the film and the skin, thus considerably reducing the friction resistance when putting on the glove. The result: The sempermed® supreme plus facilitates easy donning of the gloves.



Skin-friendly

The specially developed, automatic finishing process contains intensive and multiple washing steps. Skin-irritating substances resulting from the production process and the natural latex protein concentration are thus reduced to the technically lowest level.



Secure grip

This surgical glove offers enhanced grip for better instrument handling, thanks to its micro-rough surface and offline chlorination process.

Technical data

Type	sterile surgical glove for single use, powder-free, with synthetic lining		
Color	natural white		
Size / Overall length as per EN 455-2	5½, 6 and 6½ / 7, 7½ and 8 / 8½ and 9	6 and 6½ / 7 and 7½ / 8 and 8½	median 270 mm median 280 mm median 285 mm
Wall thickness measured in single layer in the palm area	0.19 - 0.23 mm		
Barrier performance	AQL 0.65		
Material	natural rubber latex (NR)		
Glove shape	fully anatomical with rolled rim		
Force at break as per EN 455-2	≥ 12 N, exceeding EN 455-2 (≥ 12 N instead of 9 N)		
Shelf life in original packaging if stored according to product specification	3 years		
Sterilization	radiation STERILE R		
Labeling & conformity to standards	EN 455-1/-2/-3/-4, EN ISO 21420, EN ISO 374-1 (Type B), EN ISO 374-2, EN 16523-1, EN ISO 374-4, ISO 374-5, EN 421, EN 556, ISO 2859-1, ISO 15223-1, ASTM F1671/F1671M		
Purpose	Single-use medical device class IIa according to MD Directive 93/42/EEC Single-use protective glove – PPE category III (protection against chemical substances for limited time)* according to PPE Regulation (EU) 2016/425		



Packaging

- Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in peel pack.
- In dispenser carton with sterilization indicator: 50 pairs
- In transport carton with sterilization indicator: 300 pairs

Article numbers

size 5½	822851521
size 6	822851601
size 6½	822851621
size 7	822851701
size 7½	822851721
size 8	822851801
size 8½	822851821
size 9	822851901